



Gynecology/Urology Procedure Chair

SONESTA 6202



Stille-Sonesta's 6202 procedure chair... the ultimate in positioning and comfort for you and your patient.



STILLE



Gynecology/Urology Procedure Chair

SONESTA 6202



The Sonesta 6202... Great Positioning For:

- Urodynamics
- General
- Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures

Engineered To Be Flexible:

- Locking mechanism for precise positioning
- Three motors for height, seat and back adjustments
- Easy positioning for elderly or physically challenged patients
- Wide range of adjustment, from a low of 16" (40cm) to 39.4" (100cm) high
- Engineered ergonomics. Infinitely adjustable (see specifications) for physician and patient comfort.
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable, able to support up to 375 lbs - 170kg for years of service in the most demanding practices.

The Sonesta 6202 chair built and engineered to be the ultimate in positioning for your practice.

The Sonesta 6202 is the ultimate in positioning flexibility. Packaged with our most popular accessories, the 6202 is fully equipped for a myriad of exam and treatment procedures. Three powerful (yet quiet) motors allow fast & precise positioning with minimal effort.

Standard Accessories:

3 motor control/foot plate	524-6359
1 cushion with 170mm cutout	524-6366-3
Paper roll holder	524-6326
2 fold up footrests	524-6331
2 heavy duty Rissler leg rests	524-6336
Drainage Receptacle	6325-1
2 hand grips/arm rests	524-6339-2
2 leg rest holders	524-6338
1 neck pillow	524-6346
1 full seat (w/o cut out)	524-6365

SPECIFICATIONS:

LENGTH WITH EXTENSION CUSHION	190cm	75"
LENGTH WITHOUT	143cm	56"
WIDTH	62cm	25"
HEIGHT (MIN)	40cm	16"
HEIGHT (MAX)	95cm	37"
LIFT CAPACITY	170kg	375lbs
WHEEL DIAMETER (LOCKING)	10cm	4"
NET WEIGHT	<125kg	<276lbs

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601

All Products carry the CE mark.

DISTRIBUTED BY:



Life-Tech, Inc.
www.life-tech.com
800-231-9841 / 281-491-6600

CONTACT EUROPE

Stille Surgical AB
Sundbybergsvägen 1A
SE-171 73 Solna, Sweden
PHONE: +46 8 588 580 00
FAX: +46 8 588 580 05
Email: info@stille.se
www.stille.se



Images may include some optional items. Contact your representative for available configurations. Specifications subject to change without notice.