



Fluoroscopy Procedure Table  
SONESTA

# 6210



The Stille-Sonesta 6210 fluoroscopy table... all the right moves and all the right views, for all of your C-Arm procedures.

# STILLE



# Fluoroscopy Procedure Table SONESTA 6210



## Positioned to be the best.

- Three motors offer unparalleled positioning/access options
- Radiolucent carbon fiber tabletop for head-to-toe imaging
- 4 Double wheel rolling casters offer unencumbered manoeuvrability

## Standard Accessories:

Hand control	524-6356
2 footrests	524- 6331
1 fixed patient arm rest	6339-X
1 plastic funnel w/drainage	6325-9
Detachable arm rests	524-6021-6
Armrest holders	524-6338-4
1 detachable side rail	524-6328



## The Sonesta 6210... engineered to change your point of view... simply and easily.

The 6210 fluoroscopy procedure table offers the physician standing, seated, and supine positioning capabilities. This versatile table easily adjusts from/to standing, seated, and supine positions via a hand held control... thus creating a table that is easily accessible to patients of any age. The side-mounted cantilevered design offers the physician anterior/posterior, lateral and oblique views without obstruction. The superb c-arm access and a wide range of positioning possibilities, make the 6210 the ideal table for any fluoroscopy procedure.



### SPECIFICATIONS

Length without Foot rest/leg extension:	156cm	63"
Length with Foot rest/leg extension:	196cm	78"
Width (total):	99cm	39"
Height (Min):	50cm	20"
Height (Max):	130cm	51"
Imaging area without Foot rest/leg extension (length):	1560cm	63"
Imaging area with Foot rest/leg extension (length):	196cm	78"
Weight capacity:	250kg	551lbs

### ELECTRICAL:

120/230v 24v DC Output

Meets EMC (Electromagnetic Compatibility) for IEC 60601

### DISTRIBUTED BY:



**Life-Tech<sup>®</sup>, Inc.**  
www.life-tech.com  
800-231-9841 / 281-491-6600

### CONTACT EUROPE

Stille Surgical AB  
Sundbybergsvägen 1A  
SE-171 73 Solna, Sweden  
PHONE: +46 8 588 580 00  
FAX: +46 8 588 580 05  
Email: info@stille.se  
www.stille.se

