



Gynecology/Urology Procedure Chair

SONESTA 6302



The Stille-Sonesta 6302 chair plus... dual-motor affordability and versatility for just about any type of Gynecology and/or Urology procedure.



STILLE



Gynecology/Urology Procedure Chair

SONESTA 6302

The 6302 Fits Your Needs.

- Urodynamics
- General
- Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures



Engineered To Be Flexible.

- Sonesta flexibility for office, clinic, or hospital
- Easy positioning for patient and physician comfort
- Accepts the full range of Sonesta Accessories
- Wide range of adjustment, from a low of 16" (40cm) to 39.4" (100cm) high
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable supports up to 375 lbs - 170 kg for years of service in the most demanding practices.

Standard Accessories:

- 2 motor control/foot plate
524- 6353/524- P6409
- Cushion w/170mm cutout 524-6366-3
- Paper roll holder 524-6326
- 2 heavy duty Rissler leg rests 524-6336
- Drainage Receptacle 524-6325-1
- 2 hand grips/arm rests 524-6339-2
- 2 leg rest holders 524-6338

The Sonesta 6302 engineered to be the best, most flexible solution in your examination room.

This two-motor table offers complete flexibility for offices on a budget. Motorized height & seat position, with gas piston assisted control of back angle. Or, choose the 6302U: fixed seat with three tilt-position settings and two-motored powered height and back control. Both models include a basic accessory package to cover many procedures, and they accept the full line of Sonesta accessories.

SPECIFICATIONS:

LENGTH WITH EXTENSION CUSHION	175cm	69"
LENGTH WITHOUT	125cm	49.2"
WIDTH	65cm	25.5"
HEIGHT (MIN)	40cm	16"
HEIGHT (MAX)	100cm	39"
LIFT CAPACITY	170kg	375lbs
WHEEL DIAMETER (LOCKING)	10cm	4"
NET WEIGHT	<100kg	<220lbs

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.

DISTRIBUTED BY:



CONTACT EUROPE

Stille Surgical AB
Sundbybergsvägen 1A
SE-171 73 Solna, Sweden
PHONE: +46 8 588 580 00
FAX: +46 8 588 580 05
Email: info@stille.se
www.stille.se

