



Gynecology/Urology Procedure Chair

SONESTA **6303**



The 6303...
three motor
positioning makes
this a versatile,
comfortable solution
for any Gynecology
and/or Urology
procedure.



STILLE



Gynecology/Urology Procedure Chair

SONESTA 6303

The 6303 Fits Your Needs.

- Urodynamics
- General
- Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures



Engineered To Be Flexible.

- Three Motors for height, seat and back adjustments
- Easy positioning for patient and physician comfort
- Accepts the full range of Stille Accessories
- Wide range of adjustment, from a low of 16" (40cm) to 39.4" (100cm) high
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable supports up to 375 lbs - 170 kg for years of service in the most demanding practices.

The 6303 is engineered to be the best, most flexible solution in your examination room.

The 6303 is an affordable path to Stille flexibility. Three motors provide full height, seat, and back adjustments. Packaged with basic accessories, the 6303 accepts the full line of procedure enhancing Stille add-ons, allowing you to create the perfect table for your needs. The 6303's ergonomic design delivers a higher level of patient and physician comfort and flexibility... making it an incredible value in today's market!

Standard Accessories:

- 3 motor foot control 524-6357
- Cushion w/170mm cutout 524-6366-3
- Paper roll holder 524-6326
- 2 heavy duty Rissler leg rests 524-6336
- 2 leg rest holders 524-6338
- Drainage Receptacle 524-6325-1
- 2 hand grips/arm rests 524-6339-2

SPECIFICATIONS:

LENGTH WITH EXTENSION CUSHION	175cm	69"
LENGTH WITHOUT	125cm	49.2"
WIDTH	65cm	25.5"
HEIGHT (MIN)	40cm	16"
HEIGHT (MAX)	100cm	39"
LIFT CAPACITY	170kg	375lbs
WHEEL DIAMETER (LOCKING)	10cm	4"
NET WEIGHT	<100kg	<220lbs

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.

CONTACT USA:



Life-Tech[®], Inc.
www.life-tech.com
800-231-9841 / 281-491-6600

STILLE PRODUCTS:

- Urology, Gynecology and Ultrasound Procedure Chairs
- Radiolucent-Cardio-Vascular Surgical Tables
- Surgical Instruments
- Surgical Retractor Systems

