

THE POWER OF **NOW**



NEXT GEN **DESIGN**
MEETS CLINICAL **APPLICATION**

Elle 450

- Genera has taken ideas and suggestions from physicians like you to provide the platform exquisitely made for a multitude of procedures. From general physical exams, to cystoscopy and urodynamics, Elle 450 is the right choice for your busy practice.
- Quick, easy positioning will allow procedures to flow smoothly and the contoured cushions provide relaxed patient comfort.
- In addition to traditional procedures, the Elle 450 design takes into consideration the advances in surgical techniques and is well suited, with the appropriate accessories, to bring minor operative procedures into the clinic or office.
- Make Elle 450 your choice for superior performance and value.



 **genera**
nextgenerationpositioning



SPECIFICATIONS

Height (max)	37"
Height (min)	20"
Length (w/o extension)	51"
Length (w extension)	72"
Length (back)	32"
Length (seat)	18"
Length (extension)	20"
Width (back)	23"
Width (seat)	23"
Width (extension)	21"
Table Weight	300lbs.
Patient Weight Capacity	500lbs.
Power Input	120/230V ~50/60Hz
Power Output	24VDC (to motor)

ELLE 450

Designed with motorized height and back adjustments.

(1) exam seat cushion	Included
(1) commode seat cushion	Included
(2) contoured arm rests	Included
(1) paper roll holder	Included
(2) adjustable foot rests	Included
(2) accessory rails (back)	Included
(2) accessory rails (seat)	Included

OPTIONAL ACCESSORIES

Knee crutches
Table extension
Urodynamic funnel
Others available as required



Distributed by:



Life-Tech, Inc
Life-Tech International, Inc

LIFE-TECH, INC.
4235 Greenbriar Drive
Stafford, TX 77477
Phone: 281-491-6600 Fax: 281-491-7197
www.life-tech.com

Tables are manufactured in an ISO 13485:2003 certified facility and meet and comply with ETL or UL Standards: CAN/CSA C22.2 Bi 601.1-M90-Medical Electrical Equipment-Part 1: General Requirements for Safety, including: C22.2 No. 01.1S1-94 (IEC601-1, Amendment 1:1991) Supplement No. 1-94 to CAN/CSA 22.2 No. 601.12-Mio. Also includes UL 2601-1, Second Addition (1997) "Medical Electrical Equipment, Part 1: "General Requirements for Safety" 2nd Edition including Amendments 1 and 2.



GENERA TABLES, LLC
8215 SANTA CLARA DRIVE
DALLAS, TEXAS 75218
877.321.1819 FAX: 877.826.1100
GENTERAMED.COM