

## FOLDaCUP™ COLLAPSIBLE MEASURING CUP

- Easily mailed in a standard envelope
- Calibrated accuracy<sup>(1)</sup>
- Ample capacity<sup>(2)</sup>
- Improves patient satisfaction and compliance<sup>(3)</sup>

The FOLDaCUP is designed for home use by patients for recording intake or bladder diaries. It solves the problem of getting a measuring vessel to the patient because it can be mailed or given to the patient in the same envelope as the diary form.

The FOLDaCUP is very easy to assemble and use, and its portability facilitates use away from home. It is supplied folded and packaged in a specially designed envelope that can be quickly inserted into a standard mailing envelope.

### SPECIFICATIONS

**Size:** Folded – 7”x 7”  
Assembled – 5 3/8” high X 3 3/4” wide  
X 3 1/4” deep

**Capacity:** 1,000cc

**Accuracy:** (95% confidence limits)<sup>(1)</sup>:

| 100cc | 500cc | 1,000cc |
|-------|-------|---------|
| ±9%   | ±4%   | ±3%     |
| ±9cc  | ±21cc | ±32cc   |

- (1) Based on repeat measurements of 10 randomly selected FOLDaCUPS. Data available from Life-Tech on request.
- (2) More than 17% of diaries from asymptomatic subjects record maximum voided volumes of over 700cc. Data available from Life-Tech on request.
- (3) If asked to purchase a measuring vessel, the patient will often acquire a kitchen measuring cup with insufficient capacity and variable accuracy, or may simply not purchase a cup and falsify the volume measurements.



Patent Pending



**Life-Tech®**  
Life-Tech International

P.O. Box 1849  
Stafford, TX 77497-1849  
281-491-6600 800-231-9841  
FAX: 281-491-7197  
email: [itsales@life-tech.com](mailto:itsales@life-tech.com)  
website: [www.life-tech.com](http://www.life-tech.com)

